

# ERP27 Polyclonal Antibody



Catalog Number: 162921

## Product name

ERP27 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to ERP27

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-273 of human ERP27 (NP\_689534.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



## Western blot - ERP27 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ERP27 antibody at 1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.