

Anti-CNGA1 antibody



Catalog Number: 175921

Product name

Polyclonal

Anti-CNGA1 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat

Applications

Antibody description

WB

Rabbit polyclonal antibody to CNGA1

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cnga1. aa range:401-450

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1.; Western Blot analysis of HepG2 cells using CNG-1 Polyclonal Antibody