

FDXR Polyclonal Antibody



Catalog Number: 164221

Product name

FDXR Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to FDXR

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 252-491 of human FDXR (NP_077728.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:100 - 1:200

Validations



Western blot - FDXR Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using FDXR antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 30s.


Immunohistochemistry - FDXR Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using FDXR antibody at dilution of 1:100 (40x lens).

FDXR Polyclonal Antibody




Catalog Number: 164221



Immunohistochemistry - FDXR Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophagus using FDXR antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - FDXR Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using FDXR antibody at dilution of 1:100 (40x lens).