

Anti-UBA52 Rabbit antibody



Catalog Number: 167321

Product name

Anti-UBA52 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

UBA52 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human UBA52

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000-1/5000


IHC: 1/100-1/500

ICC/IF: 1/20-1/100


FC: 1/500

IP: 1/20

Validations



Immunocytochemistry of UBA52(green) in HeLa cells using UBA52 Rabbit pAb at dilution 1/50, and



DAPI(blue)

Western blot detection of UBA52 in HeLa,SKOV-3 cell lysates using UBA52 Rabbit pAb(1:1000 diluted).Predicted band size:15KDa.Observed band size:10KDa.