

Napsin A antibody



Catalog Number: 107321

Product name

Napsin A antibody

IF: 1:10-1:100

Specificity

Human, Mouse; other species not tested.

Antibody description

Napsin A Mouse Monoclonal antibody. Positive IHC detected in human lung cancer tissue, human kidney tissue, human lung tissue, rat lung tissue. Positive IF detected in HUVEC cells. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of Napsin A recombinant protein (Accession Number: NM_004851). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, IHC, IF

Dilutions


Recommended Dilution:

IHC: 1:20-1:200


Validations



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:107321(NAPSA antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:107321(NAPSA antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HUVEC cells using Catalog No:107321(NAPSA Antibody) at dilution

Napsin A antibody



Catalog Number: 107321

of 1:25 and Rhodamine-Goat anti-Mouse IgG