

SNRNP70 Polyclonal Antibody



Catalog Number: 164721

Product name

SNRNP70 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to SNRNP70

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human SNRNP70 (NP_003080.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:1000


IHC 1:50 - 1:100

Validations




Western blot - SNRNP70 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SNRNP70 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunohistochemistry - SNRNP70 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat liver using SNRNP70 antibody at dilution of 1:100



(40x lens).


SNRNP70 Polyclonal Antibody



Catalog Number: 164721

Immunohistochemistry - SNRNP70 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat testis using SNRNP70 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - SNRNP70 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat ovary using SNRNP70 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - SNRNP70 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using SNRNP70 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - SNRNP70 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat spleen using SNRNP70 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - SNRNP70 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using SNRNP70 antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - SNRNP70 Polyclonal Antibody

SNRNP70 Polyclonal Antibody



Catalog Number: 164721


Immunohistochemistry of paraffin-embedded rat pancreas using SNRNP70 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - SNRNP70 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human tonsil using SNRNP70 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - SNRNP70 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using SNRNP70 antibody at dilution




of 1:100 (40x lens).

Immunohistochemistry - SNRNP70 Polyclonal


Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using SNRNP70 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - SNRNP70 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using SNRNP70 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - SNRNP70 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using SNRNP70 antibody at dilution of 1:100 (40x lens).

SNRNP70 Polyclonal Antibody




Catalog Number: 164721



Immunohistochemistry - SNRNP70 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human uterus using SNRNP70 antibody at dilution



of 1:100 (40x lens).


Immunohistochemistry - SNRNP70 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human uterine cancer using SNRNP70 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - SNRNP70 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse testis using SNRNP70 antibody at dilution



of 1:100 (40x lens).


Immunohistochemistry - SNRNP70 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse spinal cord using SNRNP70 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - SNRNP70 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using SNRNP70 antibody at dilution



of 1:100 (40x lens).


SNRNP70 Polyclonal Antibody



Catalog Number: 164721

Immunohistochemistry - SNRNP70 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse spleen using SNRNP70 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - SNRNP70 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using SNRNP70 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - SNRNP70 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse pancreas using SNRNP70 antibody at dilution of 1:100 (40x lens).