

SETDB2 antibody



Catalog Number: 115131

Product name

SETDB2 antibody

Specificity

Human; other species not tested.

Antibody description

SETDB2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 101 kDa

Preparation

This antibody was obtained by immunization of SETDB2 recombinant protein (Accession Number: NM_001160308). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115131(SETDB2 Antibody) at dilution of 1:300