

Anti-Phospho-PAK1(Ser144)/2(Ser141)/3 (Ser154) Rabbit antibody



Catalog Number: 167231

Product name

Anti-Phospho-PAK1(Ser144)/2(Ser141)/3 (Ser154)
Rabbit antibody

WB: 1/2000-1/10000

IHC: 1/20-1/100

ICC/IF: 1/50

FC: 1/20

Specificity

Human, Mouse, Rat

IP: 1/20

Antibody description

Phospho-PAK1(Ser144)/2(Ser141)/3 (Ser154)
Rabbit polyclonal antibody

Preparation

Antigen: A synthetic phosphopeptide
corresponding to residues surrounding Ser144 of
human PAK1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M
NaCl, 40%Glycerol, 0.01% sodium azide and
0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of
receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions

Validations



Western blot detection of PAK1 (Phospho-Ser144)/PAK2 (Phospho-Ser141)/PAK3 (Phospho-Ser154) in C6,3T3,HeLa cell lysates using PAK1 (Phospho-Ser144)/PAK2 (Phospho-Ser141)/PAK3 (Phospho-Ser154) Rabbit pAb(1:1000 diluted).Predicted band size:61kDa.Observed band size:61kDa.



Immunocytochemistry of PAK1 (Phospho-Ser144)/PAK2 (Phospho-Ser141)/PAK3 (Phospho-Ser154)(green) in HeLa cells using PAK1 (Phospho-Ser144)/PAK2 (Phospho-Ser141)/PAK3 (Phospho-Ser154) Rabbit pAb at dilution 1/50, and DAPI(blue)