

# Anti-Hsp27 (7E5) Mouse antibody



Catalog Number: 168231

## Product name

Anti-Hsp27 (7E5) Mouse antibody

## Specificity

Human, Monkey

## Antibody description

Hsp27 (7E5) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human Hsp27 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, ICC/IF

## Dilutions

WB: 1/1000


;ICC: 1/300

## Validations



Western blot detection of Hsp27 in HeLa, HsHela, A431, COS7, Lncap and MCF7 cell lysates and using Hsp27 mouse mAb (1:1000 diluted). Predicted band size: 23KDa. Observed

band size: 27KDa.



Immunocytochemistry stain of HeLa using Hsp27 mouse mAb (1:300).