

# Anti-STAT5a/b Rabbit antibody



Catalog Number: 168631

## Product name

Anti-STAT5a/b Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

STAT5a/b Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human STAT5A

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/1000

ICC/IF: 1/50


FC: 1/20

IP: 1/50

## Validations



Immunocytochemistry of STAT5A/B(green) in HeLa cells using STAT5A/B Rabbit mAb at dilution 1/50, and DAPI(blue)



and DAPI(blue)

Western blot detection of STAT5a/b in K562,C6,3T3,HeLa cell lysates using STAT5a/b Rabbit mAb(1:1000 diluted).Predicted band size:90kDa.Observed band size:90kDa.