

# Anti-Chk1 (6C5) Mouse antibody



Catalog Number: 168441

## Product name

Anti-Chk1 (6C5) Mouse antibody

## Specificity

Human, Mouse

## Antibody description

Chk1 (6C5) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of CHK1 expressed in E. Coli.

## Formulation

Ascitic fluid containing 0.03% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB

## Dilutions


WB: 1/500 - 1/2000

ELISA: 1/10000

## Validations



Western blot analysis of extracts from HeLa cells, untreated or treated with UV, using Chk1 (Phospho-Ser280) Rabbit pAb (310037, 1:500 diluted, upper) or CHK1 Mouse mAb



(168441, 1:500 diluted, lower).

Western blot analysis using CHK1 mouse mAb against A431 (1), HeLa (2), NIH/3T3 (3) and K562 (4) cell lysate.