

# Timp-3 antibody

Catalog Number: 116141

---

**Product name**

Timp-3 antibody

**Specificity**

Human, Mouse, Rat; other species not tested.

**Antibody description**

Timp-3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human kidney tissue, human lung cancer tissue, human oesophagus cancer tissue, human stomach tissue.

**Preparation**

This antibody was obtained by immunization of Timp-3 recombinant protein (Accession Number: NM\_000362). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, IHC


**Dilutions**

Recommended Dilution:

IHC: 1:20-1:200

**Validations**

Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116141(Timp-3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116141(Timp-3 Antibody) at dilution of 1:50 (under 40x lens)