

ARD1B antibody



Catalog Number: 108241

Product name

Polyclonal

ARD1B antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

ARD1B Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 26kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of ARD1B recombinant protein (Accession Number: BC004552). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108241(ARD1B antibody) at dilution of 1:500