

# Anti-FDFT1 Rabbit antibody

Catalog Number: 168641

## Product name

Anti-FDFT1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

FDFT1 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human FDFT1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, IP


## Dilutions

WB: 1/1000


ICC/IF: 1/50

IP: 1/20

## Validations



Western blot detection of FDFT1 in K562, C6, HeLa cell lysates using FDFT1 Rabbit mAb(1:1000 diluted). Predicted band size: 48kDa. Observed band size: 48kDa.



Immunocytochemistry of FDFT1(green) in HeLa cells using FDFT1 Rabbit mAb at dilution 1/200, and DAPI(blue)