

CAMK1G Polyclonal Antibody



Catalog Number: 165341

Product name

CAMK1G Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CAMK1G

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 217-476 of human CAMK1G (NP_065172.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - CAMK1G Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CAMK1G antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.




Western blot - CAMK1G Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CAMK1G antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

CAMK1G Polyclonal Antibody




Catalog Number: 165341



Immunohistochemistry - CAMK1G Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using CAMK1G Antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - CAMK1G Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human brain cancer using CAMK1G Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - CAMK1G Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using CAMK1G Antibody at dilution of 1:100 (40x lens).