

POLR3F antibody



Catalog Number: 114051

Product name

POLR3F antibody

Immunogen

[Human POLR3F Recombinant protein \(GST tag & His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

POLR3F Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, BxPC-3 cells, HeLa cells, mouse skin tissue. Observed molecular weight by Western-blot: 36-39kd

Preparation

This antibody was obtained by immunization of POLR3F recombinant protein (Accession Number: NM_006466). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114051(POLR3F antibody) at dilution of 1:400