

Anti-IGLC2 antibody



Catalog Number: 175051

Product name

Anti-IGLC2 antibody

Specificity

Human

Antibody description

Mouse monoclonal antibody to IGLC2

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*TBS (pH7.4), 1%BSA, 40%Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P

Dilutions

WB: 1:500

IHC-P: 1:50-1:200

Validations






Fig1: Western blot analysis of IGLC2 on human IGLC2 recombinant protein using anti-IGLC2



antibody at 1/1,000 dilution.

Fig2: Western blot analysis of IGLC2 on HEK293 (1) and IGLC2-hIgGFc transfected HEK293 (2) cell lysate using anti-IGLC2 antibody at 1/1,000



dilution.

Fig3: Immunohistochemical analysis of paraffin-embedded human colon tissues using anti-IGLC2 antibody. Counter stained with hematoxylin.

Anti-IGLC2 antibody

Catalog Number: 175051




Fig4: Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissue using anti-IGLC2 antibody. Counter stained with hematoxylin.