

SAMD14 antibody



Catalog Number: 115051

Product name

SAMD14 antibody

Specificity

Human; other species not tested.

Antibody description

SAMD14 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human colon cancer tissue. Positive WB detected in COLO 320 cells, HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 50-55 kDa

Preparation

This antibody was obtained by immunization of SAMD14 recombinant protein (Accession Number: NM_001257359). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:115051(SAMD14 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:115051(SAMD14 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:115051(SAMD14 Antibody) at dilution of 1:50

SAMD14 antibody



Catalog Number: 115051

(under 40x lens)