

CD200 antibody



Catalog Number: 107051

Product name

CD200 antibody

IHC: 1:50-1:500


Specificity

Human; other species not tested.

Validations

Antibody description

CD200 Mouse Monoclonal antibody. Positive WB detected in human spleen tissue, fetal human brain tissue. Positive IHC detected in human tonsillitis tissue. Observed molecular weight by Western-blot: 45-50 kDa




Preparation

This antibody was obtained by immunization of CD200 recombinant protein (Accession Number: NM_005944). Purification method: Protein G purified.

human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:107051(CD200 Antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Monoclonal

Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using (CD200 Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).

Ig Type

Mouse IgG1



Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using (CD200

CD200 antibody



Catalog Number: 107051

Antibody) at dilution of 1:200 (under 40x lens).

heat mediated antigen retrieved with Tris-EDTA
buffer(pH9).