

CD134/OX40 antibody



Catalog Number: 109051

Product name

CD134/OX40 antibody

Specificity

Human; other species not tested.

Antibody description

CD134/OX40 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, K-562 cells, L02 cells, U-937 cells. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003327).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109051(TNFRSF4 antibody) at dilution of 1:300