

HSD17B8 antibody



Catalog Number: 111551

Product name

HSD17B8 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HSD17B8 Rabbit Polyclonal antibody. Positive IHC detected in human liver. Positive WB detected in human testis tissue, mouse small intestine tissue. Observed molecular weight by Western-blot: 34kd

Preparation

This antibody was obtained by immunization of HSD17B8 recombinant protein (Accession Number: NM_014234). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:111551(HSD17B8 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver using Catalog No:111551(HSD17B8 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:111551(HSD17B8 antibody) at dilution of 1:100 (under 40x lens)