

# UBL4A antibody



Catalog Number: 116651

## Product name

UBL4A antibody

WB: 1:200-1:1000

IP: 1:200-1:2000

## Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

IF: 1:20-1:200

## Antibody description

UBL4A Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive IF detected in HepG2 cells. Positive WB detected in HeLa cells, human brain tissue, mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 17 kDa

## Preparation

This antibody was obtained by immunization of UBL4A recombinant protein (Accession Number: NM\_014235). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IF, IP


## Dilutions

Recommended Dilution:

## Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116651(UBL4A antibody) at dilution of 1:400




Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:116651(UBL4A antibody) at dilution of 1:100 (under 10x lens)

# UBL4A antibody


Catalog Number: 116651



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:116651(UBL4A antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using UBL4A antibody Catalog No:116651 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent



DNA dye).

IP Result of anti-UBL4A (IP:Catalog No:116651, 4ug; Detection:Catalog No:116651 1:500) with mouse testis tissue lysate 4400ug.