

Anti-MEK3 Rabbit antibody



Catalog Number: 168651

Product name

Anti-MEK3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

MEK3 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human MEK3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/50


FC: 1/100

IP: 1/20

Validations



Western blot detection of MEK3 in C6, 3T3, HeLa cell lysates using MEK3 Rabbit mAb (1:1000 diluted). Predicted band size: 39 kDa. Observed band size: 38-40 kDa.



Immunocytochemistry of MEK3 (green) in HeLa cells using MEK3 Rabbit mAb at dilution 1/50, and DAPI (blue)