

PXR antibody



Catalog Number: 114351

Product name

ELISA, WB, IP

PXR antibody

Dilutions

Specificity

Recommended Dilution:

Human, Mouse, Rat; other species not tested.

WB: 1:200-1:2000

IP: 1:200-1:2000


Antibody description

PXR Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, A549 cells, MCF7 cells, mouse colon tissue. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 45-50kd

Validations

Preparation

This antibody was obtained by immunization of PXR recombinant protein (Accession Number: NM_003889). Purification method: Antigen affinity purified.



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114351(NR1I2 antibody) at dilution of 1:200

Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications



IP Result of anti-PXR (IP:Catalog No:114351, 4ug; Detection:Catalog No:114351 1:500) with HepG2 cells lysate 3600ug.