

# ATP5C1 antibody



Catalog Number: 108351

## Product name

ATP5C1 antibody

## Clonality

Polyclonal

## Immunogen

[Human ATP5C1 Recombinant protein \(GST tag & His tag\)](#)

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

ATP5C1 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, human heart tissue, mouse brain tissue. Observed molecular weight by Western-blot: 33 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of ATP5C1 recombinant protein (Accession Number: XM\_047425254). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:108351(ATP5C1 antibody) at dilution of 1:500