

# DPCD antibody



Catalog Number: 110061

## Product name

DPCD antibody

## Immunogen

[Human RP11-529I10.4 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DPCD Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue. Positive WB detected in human testis tissue, mouse testis tissue. Observed molecular weight by Western-blot: 23-26 kDa

## Preparation

This antibody was obtained by immunization of DPCD recombinant protein (Accession Number: NM\_015448). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110061(RP11-529I10.4 antibody) at dilution of



1:300

Immunohistochemistry of paraffin-embedded human heart slide using Catalog No:110061(RP11-529I10.4 Antibody) at dilution of 1:50