

Anti-Bcl10 Rabbit antibody



Catalog Number: 168061

Product name

Anti-Bcl10 Rabbit antibody

Specificity

Human

Antibody description

Bcl10 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Bcl10

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

IP: 1/20

Validations



Western blot detection of Bcl10 in A549,HL-60,U2OS cell lysates using Bcl10 Rabbit mAb(1:1000 diluted).Predicted band size:26kDa.Observed band size:26kDa.