

FAM116B antibody



Catalog Number: 110461

Product name

FAM116B antibody

Immunogen

[Human FAM116B Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

FAM116B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A549 cells. Observed molecular weight by Western-blot: 66 kDa

Preparation

This antibody was obtained by immunization of FAM116B recombinant protein (Accession Number: NM_001001794). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog

No:110461(FAM116B Antibody) at dilution of 1:600