

GAN antibody



Catalog Number: 110861

Product name

Polyclonal

GAN antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

GAN Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HEK-293 cells. Observed molecular weight by Western-blot: 68kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of GAN recombinant protein (Accession Number: NM_022041). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110861(GAN antibody) at dilution of 1:500