

ACADM Polyclonal Antibody



Catalog Number: 161861

Product name

Rabbit IgG

ACADM Polyclonal Antibody

Applications**Specificity**

WB, IHC

Human, Mouse

Dilutions**Antibody description**


WB 1:500 - 1:2000

Polyclonal antibody to ACADM

IHC 1:50 - 1:100

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 26-421 of human ACADM (NP_000007.1).

Validations**Formulation**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Immunohistochemistry - ACADM Polyclonal Antibody

Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human kidney using ACADM Antibody at dilution of 1:200 (40x lens).

Ig Type