

TNS4 antibody



Catalog Number: 116161

Product name

TNS4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TNS4 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Positive FC detected in HepG2 cells. Positive IHC detected in human liver cancer tissue. Observed molecular weight by Western-blot: 77kd

Preparation

This antibody was obtained by immunization of TNS4 recombinant protein (Accession Number: BC018706). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, FC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:116161(TNS4 antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:116161(TNS4 Antibody) at dilution of 1:50 (under 40x lens)



1X10⁶ HepG2 cells were stained with .2ug TNS4 antibody (Catalog No:116161, red) and control antibody (blue). Fixed with 90% MeOH blocked

TNS4 antibody



Catalog Number: 116161

with 3% BSA (30 min). Alexa Fluor 488-

conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)
with dilution 1:1000.