

Fibulin 5 antibody



Catalog Number: 110661

Product name

Fibulin 5 antibody

WB: 1:200-1:2000

IP: 1:200-1:1000

Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

Antibody description

Fibulin 5 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue, human heart tissue, human pancreas tissue. Positive IP detected in mouse heart tissue. Positive WB detected in mouse heart tissue, mouse placenta tissue, rat brain tissue. Observed molecular weight by Western-blot: 64-66 kDa

Validations



Preparation

This antibody was obtained by immunization of Fibulin 5 recombinant protein (Accession Number: NM_006329). Purification method: Antigen affinity purified.

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110661(FBLN5 antibody) at dilution of 1:400

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:



IP Result of anti-FBLN5 (IP:Catalog No:110661, 3ug; Detection:Catalog No:110661 1:300) with mouse heart tissue lysate 8000ug.



Fibulin 5 antibody



Catalog Number: 110661

Immunohistochemical of paraffin-embedded
human lung cancer using Catalog

No:110661(FBLN5 antibody) at dilution of 1:50
(under 10x lens)