

GAS41 antibody



Catalog Number: 110871

Product name

Polyclonal

GAS41 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

GAS41 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, A549 cells. Observed molecular weight by Western-blot: 27-30 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of GAS41 recombinant protein (Accession Number: NM_006530). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110871(YEATS4 antibody) at dilution of 1:300