

Anti-MEKK2 Rabbit antibody



Catalog Number: 166871

Product name

Anti-MEKK2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

MEKK2 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human MEKK2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/2000-1/10000

ICC/IF: 1/50


FC: 1/20

IP: 1/20

Validations



Immunocytochemistry of MEKK2(green) in HeLa cells using MEKK2 Rabbit pAb at dilution 1/50,



and DAPI(blue)

Western blot detection of MEKK2 in C6,3T3,HeLa cell lysates using MEKK2 Rabbit pAb(1:1000 diluted).Predicted band size:70kDa.Observed band size:70kDa.