

SGMS1 antibody



Catalog Number: 115171

Product name

SGMS1 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SGMS1 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive WB detected in human heart tissue. Positive IP detected in mouse heart tissue. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of SGMS1 recombinant protein (Accession Number: XM_047424974). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115171(SGMS1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human liver using Catalog No:115171(SGMS1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:115171(SGMS1 antibody) at dilution of 1:50 (under 10x lens)

SGMS1 antibody



Catalog Number: 115171

antibody) at dilution of 1:50 (under 40x lens)

IP Result of anti-SGMS1 (IP:Catalog No:115171, 4ug; Detection:Catalog No:115171 1:300) with mouse heart tissue lysate 4000ug.

