

# Anti-Galectin 3 (7G9) Mouse antibody



Catalog Number: 168171

## Product name

Anti-Galectin 3 (7G9) Mouse antibody

## Specificity

Human, Mouse

## Antibody description

Galectin 3 (7G9) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant protein expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB, IHC-P, ELISA


## Dilutions

WB: 1/1;000-2;000 IHC1/100-200 ELISA 10;000-50;000


## Validations



Immunohistochemical analysis of paraffin-embedded human-colon2 using antibody diluted at



Western blot detection of Galectin-3/LGALS3 in CHO-K1,Hela cell lysates using Galectin-3/LGALS3 (7G9) Mouse mAb(1:1000 diluted).Predicted band size:26KDa.Observed band size:26KDa.



Immunohistochemical analysis of paraffin-embedded human-colon using antibody diluted at

# Anti-Galectin 3 (7G9) Mouse antibody



Catalog Number: 168171

---

Western blot analysis of 1)MCF7, 2) 3T3 with  
Galectin-3 Mouse mAb diluted at 1:2,000.

