

# Anti-IL-1RAcP Rabbit antibody



Catalog Number: 167571

## Product name

Rabbit IgG

Anti-IL-1RAcP Rabbit antibody

## Applications

## Specificity

WB, FC, IP

Human

## Dilutions

## Antibody description

WB: 1/1000

IL-1RAcP Rabbit monoclonal antibody

FC: 1/20

## Preparation

IP: 1/20

Antigen: Recombinant protein of human IL-1RAcP

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA


## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of IL-1RAcP in HL-60,U251,U87-MG cell lysates using IL-1RAcP Rabbit mAb(1:1000 diluted).Predicted band size:65kDa.Observed band size:78kDa.