

TMOD3 antibody



Catalog Number: 116271

Product name

TMOD3 antibody

Immunogen

[Human TMOD3 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TMOD3 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human heart tissue, human kidney tissue, human lung tissue, mouse brain tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of TMOD3 recombinant protein (Accession Number: XM_017022088). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116271(TMOD3 Antibody) at dilution of 1:1000