

Anti-TXNRD1 Rabbit antibody



Catalog Number: 167271

Product name

Anti-TXNRD1 Rabbit antibody

FC: 1/20-1/100

IP: 1/20

Specificity

Human, Mouse, Rat

Validations

Antibody description

TXNRD1 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human TXNRD1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions


WB: 1/1000-1/5000

IHC: 1/20


ICC/IF: 1/20



Western blot detection of TXNRD1 in K562, HeLa cell lysates using TXNRD1 Rabbit pAb(1:1000 diluted). Predicted band size: 67kDa. Observed band size: 55kDa.



Immunohistochemistry of TXNRD1 in paraffin-embedded Human lung cancer tissue using TXNRD1 Rabbit pAb at dilution 1/20



Immunocytochemistry of TXNRD1(green) in HeLa cells using TXNRD1 Rabbit pAb at dilution 1/50, and DAPI(blue)