

FKBP1A antibody



Catalog Number: 110671

Product name

FKBP1A antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FKBP1A Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human heart tissue. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 12-16 kDa

Preparation

This antibody was obtained by immunization of FKBP1A recombinant protein (Accession Number: NM_054014). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:110671(FKBP1A antibody) at dilution of 1:100



(under 25x lens)

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110671(FKBP1A Antibody) at dilution of 1:1000