

Anti-SERCA2 ATPase Rabbit antibody



Catalog Number: 167371

Product name

Anti-SERCA2 ATPase Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

SERCA2 ATPase Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human SERCA2 ATPase

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/2000-1/10000

IHC: 1/20


ICC/IF: 1/50

FC: 1/20-1/100


Validations



Immunohistochemistry of SERCA2 ATPase in paraffin-embedded Human colon cancer tissue using SERCA2 ATPase Rabbit pAb at dilution 1/50



Western blot detection of SERCA2 ATPase in K562, Hela cell lysates using SERCA2 ATPase Rabbit pAb(1:1000 diluted). Predicted band size: 115kDa. Observed band size: 115kDa.



Immunocytochemistry of SERCA2 ATPase (green) in Hela cells using SERCA2 ATPase Rabbit pAb at dilution 1/50, and DAPI (blue)