

ZNF48 antibody



Catalog Number: 117181

Product name

ZNF48 antibody

Immunogen

[Human ZNF48 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZNF48 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of ZNF48 recombinant protein (Accession Number: NM_152652). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:117181(ZNF48 Antibody) at dilution of 1:600