

Anti-Lyn (1C8) Mouse antibody



Catalog Number: 168181

Product name

Anti-Lyn (1C8) Mouse antibody

Specificity

Human

Antibody description

Lyn (1C8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Lyn protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Lyn in Ramos, Raji, K562 and HeLa cell lysates using Lyn mouse mAb(dilution 1:1000). Predicted band size: 56/53 kDa. Observed band size: 56/53 kDa.