

Anti-NDUFB9 Rabbit antibody



Catalog Number: 168581

Product name 1:1000

Anti-NDUFB9 Rabbit antibody 1:50

Specificity 1:50

Human, Mouse, Rat 1:50

Antibody description 1:20

NDUFB9 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human NDUFB9

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions

Validations



Western blot detection of NDUFB9 in Rat Brain lysates using NDUFB9 Rabbit mAb(1:1000 diluted). Predicted band size:22kDa. Observed band size:22kDa.



Immunocytochemistry of NDUFB9(green) in HeLa cells using NDUFB9 Rabbit mAb at dilution 1/200, and DAPI(blue)