

CRTC2,TORC2 antibody



Catalog Number: 109581

Product name

CRTC2,TORC2 antibody

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human,Mouse,Rat; other species not tested.

Antibody description

CRTC2,TORC2 Rabbit Polyclonal antibody.
Positive IF detected in Hela cells. Positive IHC detected in human gliomas tissue. Positive WB detected in HepG2 cells, HEK-293 cells, mouse liver tissue. Observed molecular weight by Western-blot: 73kd

Validations



Preparation

This antibody was obtained by immunization of CRTC2,TORC2 recombinant protein (Accession Number: NM_181715). Purification method: Antigen affinity purified.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109581(CRTC2,TORC2 antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF


Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:109581(CRTC2,TORC2 antibody) at dilution of 1:100 (under 10x lens)



CRTC2,TORC2 antibody



Catalog Number: 109581

Immunofluorescent analysis of Hela cells, using CRTC2 antibody Catalog No:109581 at 1:50 dilution and FITC-labeled donkey anti-rabbit IgG

(green). Blue pseudocolor = DAPI (fluorescent DNA dye).