

# Anti-CD44 Rabbit antibody



Catalog Number: 166981

## Product name

Anti-CD44 Rabbit antibody

## Specificity

Human

## Antibody description

CD44 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human CD44

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC


## Dilutions

WB: 1/1000


IHC: 1/20

FC: 1/20


## Validations



Immunohistochemistry of CD44 in paraffin-embedded Human tonsil using CD44 Rabbit pAb



Western blot detection of CD44 in HeLa cell lysates using CD44 Rabbit pAb(1:1000 diluted).Predicted band size:82kDa.Observed band



size:82kDa.

Western blot detection of CD44 in HeLa cell lysates using CD44 Rabbit pAb(1:1000 diluted).Predicted band size:82KDa.Observed band size:82KDa.