

Anti-Cullin 1/CUL-1 Rabbit antibody



Catalog Number: 168681

Product name

WB, ICC/IF, FC

Anti-Cullin 1/CUL-1 Rabbit antibody

Dilutions

Specificity

1:1000

Human, Mouse, Rat

1:20-1:50

Antibody description

1:20

Cullin 1/CUL-1 Rabbit monoclonal antibody

Validations

Preparation

Antigen: A synthetic peptide of human Cullin 1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.


Clonality

Monoclonal


Ig Type

Rabbit IgG

Applications



Immunocytochemistry of Cullin 1(green) in HeLa cells using Cullin 1 Rabbit mAb at dilution 1/200, and DAPI(blue)



and DAPI(blue)

Western blot detection of Cullin 1/CUL-1 in K562, C6 cell lysates using Cullin 1/CUL-1 Rabbit mAb(1:1000 diluted). Predicted band size: 90kDa. Observed band size: 90kDa.