

Anti-FGFR4/CD334 antibody



Catalog Number: 103381

Product name

Anti-FGFR4/CD334 antibody

Immunogen

[Human FGFR4/CD334 \(His Tag\) recombinant protein](#)

Specificity

Human FGFR4 / CD334

Antibody description

Rabbit polyclonal to FGFR4/CD334

Preparation

Produced in rabbits immunized with purified, recombinant Human FGFR4 / CD334 (rh FGFR4 / CD334; NP_002002.3; Met 1-Asp 369). FGFR4 / CD334 specific IgG was purified by Human FGFR4 / CD334 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human FGFR4. The detection limit for Human FGFR4 is approximately 0.00975 ng/well.


IHC-P: 0.1-2 µg/mL

Validations



FGFR4 / CD334 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human FGFR4 in human kidney (from 2 donors) with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).



FGFR4 / CD334 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human FGFR4 in human stomach (from 3 donors) with rabbit

Anti-FGFR4/CD334 antibody



Catalog Number: 103381

polyclonal antibody (1 µg/mL, formalin-fixed

paraffin embedded sections).