

# TMEM106B antibody



Catalog Number: 116091

## Product name

TMEM106B antibody

WB: 1:500-1:5000

IP: 1:500-1:5000

## Specificity

Human, Mouse; other species not tested.

IHC: 1:20-1:200

IF: 1:20-1:200

## Antibody description

TMEM106B Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, A549 cells, HeLa cells, human placenta tissue. Positive IP detected in mouse brain tissue. Positive IHC detected in human brain tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 31-35kd

## Preparation

This antibody was obtained by immunization of TMEM106B recombinant protein (Accession Number: NM\_001134232). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IF, IP


## Dilutions

Recommended Dilution:

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116091(TMEM106B antibody) at dilution of 1:2000




Immunohistochemical of paraffin-embedded human brain using Catalog No:116091(TMEM106B antibody) at dilution of 1:100 (under 10x lens)

# TMEM106B antibody




Catalog Number: 116091



Immunohistochemical of paraffin-embedded human brain using Catalog No:116091(TMEM106B antibody) at dilution of 1:100 (under 40x lens)

IP Result of anti-TMEM106B (IP:Catalog No:116091, 3ug; Detection:Catalog No:116091 1:1000) with mouse brain tissue lysate 6000ug.



Immunofluorescent analysis of HeLa cells using Catalog No:116091(TMEM106B Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

