

PAH antibody



Catalog Number: 113491

Product name

PAH antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PAH Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse kidney tissue. Positive IHC detected in human liver cancer tissue. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of PAH recombinant protein (Accession Number: NM_000277). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:50-1:500

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113491(PAH Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:113491(PAH Antibody) at dilution of 1:200 (under 40x lens).